
DAS ESCAPEE PHÄNOMEN

Copyright Becton-Dickinson

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

Figure Legend (1. Abbildung: CD8 von Hersteller 1, 2. Abbildung: CD8 von Hersteller 2): Region 1 includes the CD8-positive and CD8-negative single lymphocytes within the lymphocyte gate. Region 2 contains the monocyte population. Region 3 includes a population of bright CD8-positive events that are outside of the standard lymphocyte gate.

Use of various monoclonal antibodies such as CD8 and CD41 may cause cell aggregation in immunophenotyping assays, resulting in lower percentage positive values for the antigen of interest. The number of aggregates can vary

depending on whether a FITC- or PE-conjugated monoclonal antibody is used. When analyzing lymphocyte-gated CD8-positive cells, CD8 FITC conjugates have been shown to yield lower percentage positive values than CD8 PE conjugates for the same donor.¹ Escapees can also be detected when using an unconjugated purified monoclonal antibody and a Goat Anti-Mouse (GAM) Ig FITC second-step reagent. The number of escapees is partially dependent on the cell preparation method used. A peripheral blood mononuclear cell separation method will result in detection of more escapees than will a whole blood preparation stained and then lysed with FACS Lysing Solution.

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To decrease the number of escapees present, Becton Dickinson recommends using whole blood preparations that have been lysed with FACS Lysing Solution prior to flow cytometric analysis. See the FACS Lysing Solution package insert for a detailed procedure. For cell suspensions prepared by Ficoll-Paque separation, it is recommended that a high-protein medium be used, ie PBS containing 2% fetal calf serum and 0.1% sodium azide. Cell suspensions should be vortexed at a low speed, and samples should not be centrifuged with forces greater than 300 x g. Rough handling of cells such as excessive washing or high-speed vortexing may cause increased cell aggregation. Avoid the use of agents for monocyte depletion by phagocytosis (iron carbonyl particles) as these may increase the incidence of escapees. Even when taking these precautions during sample preparation, cell aggregation may still occur. For more information, contact the Customer Support Center at your local Becton Dickinson office.

LITERATUR

1. Jackson A. Basic Phenotyping of Lymphocytes: Selection and Testing of Reagents and Interpretation of Data. Clinical Immunology Newsletter. 1990;10(4):43-55.

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